

PROCESS AND DEVICE FOR DETERMINATION OF CELL VIABILITY

This application claims the benefit of United States provisional application number 60/531,848, filed December 22, 2003, the entire contents of which are incorporated herein by reference. Throughout this application various publications are referenced.

- 5 The disclosures of these publications in their entireties are hereby incorporated by reference into this application in order to more fully describe the state of the art to which this invention pertains.

TECHNICAL FIELD OF THE INVENTION

- 10 The invention relates generally to a process and device for the determination of cell viability.

BACKGROUND OF THE INVENTION

- Various applications including the use of tissue and cell culture to produce materials for treatment of various health problems creates the necessity to quantitatively assess the viability of cells and tissues at all stages of growth and propagation. Currently this is accomplished using dye exclusion and other approaches that are messy, time-consuming, expensive, difficult to automate and perform with rapid throughput, not easy to quantitate and, in the final analysis, not very precise and accurate. The process of dye exclusion can also consume or modify some of the cells and, in the case of stem cells, the loss of the cells themselves is an extreme disadvantage. Thus, there exists a need for an improved process to determine cell and tissue viability.
- 15
20

SUMMARY OF THE INVENTION

The invention provides a method to determine the viability of cells by measuring the absolute and relative rate of metabolic activity and/or integrity of the cell membrane through the use of vibrational spectroscopy.

In one embodiment, the method comprises obtaining a container loaded with deuterated materials, introducing cells into the container whereby the cells are in contact with the deuterated materials, and obtaining vibrational spectra emitted by the cells. The vibrational spectra emitted by the cells are indicative of metabolism, thereby providing an indication of viability of the cells, such that greater metabolic activity is indicative of greater viability.

The vibrational spectra can be Raman spectra, infrared and/or near infrared spectra. The deuterated material can be, for example, selected from the group consisting of α -D-glucose; 6,6- dideutero- α -D-glucose; D_2O ; 3-O-methylglucose; 6,6- dideutero- α -D-3-O-trideuteromethylglucose; 6,6- dideutero- α -D-3-O-methylglucose; 6,6- dideutero - α -D-2-O-methylglucose; and deuterated amino acids.

In some embodiments, the method further comprises normalizing the Raman spectra obtained by comparing the Raman spectra obtained at a target wavenumber to the Raman spectra obtained at a reference wavenumber. Representative target wavenumbers include 960 cm^{-1} , 1270 cm^{-1} and $2400\text{-}2600\text{ cm}^{-1}$. Typical reference wavenumbers include the amide I Raman feature ($1600\text{-}1700\text{ cm}^{-1}$), the CH_2 Raman feature (1450 cm^{-1}), the amide III Raman feature ($1200\text{-}1350\text{ cm}^{-1}$), the CH stretching region ($2900\text{-}3000\text{ cm}^{-1}$), and the integral of all the Raman features ($300\text{-}1850\text{ cm}^{-1}$). In some embodiments, the Raman spectra are normalized by comparing the Raman spectra obtained at a target wavenumber to the fluorescence generated by a Raman excitation source, such as the H_2O fluorescence (980 nm).

In another embodiment, the method of determining cell viability comprises obtaining a container loaded with deuterated materials, introducing cells into the container whereby the cells are in contact with the deuterated materials, obtaining vibrational spectra emitted by the cells, placing an aliquot of the cells into a medium free of deuterated materials, obtaining vibrational spectra emitted by the cells in the non-deuterated medium, and determining the rate of decrease of emitted vibrational spectra. The vibrational spectra are indicative of metabolism, thereby providing an indication of

viability of the cells, such that a faster rate of decrease of emitted vibrational spectra is indicative of greater viability.

BRIEF DESCRIPTION OF THE FIGURES

Figure 1 shows spectral counts as a function of Raman shift (in cm^{-1}) measured from human leukocytes grown in culture and exposed to phosphate buffered saline for 30-90 minutes. These spectra serve as a control for comparison to those shown in Figure 2.

Figure 2 shows spectral counts as a function of Raman shift (in cm^{-1}) measured from human leukocytes grown in culture and exposed to phosphate buffered saline containing D_2O instead of H_2O for 30-90 minutes.

DETAILED DESCRIPTION

One of the most basic indicators of viability is the metabolism of glucose and the active maintenance of membrane integrity. A healthy cell keeps certain materials on the inside and other materials on the outside of its membranes or cell wall (in the case of plant cells). While it is possible to synthesize vesicles and micelles and other objects that exclude, to a greater or lesser extent, certain materials from their inner volume, it is not possible, using such entities, to mimic the full range of kinetics and thermodynamic behavior of living, i.e. metabolizing, cells or tissues. Having the ability to noninvasively monitor the consumption of nutrients and/or production of metabolites or waste allows a temporally continuous, direct measure of cell viability from the moment of harvest from the primary source, e.g. stem cells from umbilical cord blood, to the final growth stages of some particular daughter culture. A highly optimized and novel approach to accomplishing this type of measurement system is provided by this invention.

Minimization of sample preparation is a worthy goal because variability and inefficiency at that stage of analysis affects accuracy, precision, throughput and the cost structure of the overall process. Given the turbidity of cultures at all stages of handling, having a method that is tolerant of turbidity is desirable. This feature, and the inherent need to

work in what is always an essentially aqueous environment, discourages analytical approaches based on absorption spectroscopy. UV-visible spectroscopy has little specificity (without the addition of what are essentially invasive fluorophores) and mid-infrared spectroscopy has solvent absorption limitations. Absorption in the near IR is
5 feasible for some analytes, but the isolation and unambiguous association of spectral features with specific analytes has proven to be challenging at millimolar concentrations for tissue-like samples.

Definitions

All scientific and technical terms used in this application have meanings commonly used
10 in the art unless otherwise specified. As used in this application, the following words or phrases have the meanings specified.

As used herein, "vibrational spectroscopy" refers to spectroscopic techniques known in the art to be based on vibrational features, including Raman, resonance Raman, infrared (IR) and near infrared (NIR) spectroscopy.

15 As used herein, "Raman spectra associated with" a given component refers to those emitted Raman spectra that one skilled in the art would attribute to that component. One can determine which Raman spectra are attributable to a given component by irradiating that component in a relatively pure form, and collecting and analyzing the Raman spectra emitted by the component in the relative absence of other components. Those skilled in
20 the art are aware of available libraries that catalog known Raman spectra.

As used herein, "metabolism" refers to the physical and chemical processes occurring within a living cell that are necessary for the maintenance of life, including anabolism and catabolism, and the processes by which a particular substance is handled (as by assimilation and incorporation or by detoxification and excretion).

25 As used herein, unless context clearly indicates otherwise, "a" means at least one, and can include a plurality.

Determination of Cell Viability

Viability can be assessed using spectroscopy in many ways. One example to is to use standard chemometric techniques to integrate the observed Raman activity in and around the vicinity of 2200-2400 cm^{-1} Raman shift (Stokes) (the C-D stretching modes) and the activity in and around 1600-1700 cm^{-1} (amide I). The ratio of these two integrated signals is zero before the cells have absorbed any of the "deuterated viability agent" and increases when the cells are incubated with the deuterated agent. The ratio has a different time course depending on the specific nature of the agent used and the viability of the cells. Utilizing glucose as the agent, for example, a population of more viable cells takes up the agent faster than a population of less viable cells and slower than a population of cancerous cells.

The invention provides a method to determine the viability of cells by measuring the absolute and relative rate of metabolic processes. In one embodiment, the method measures indicia of glucose uptake, clearing and/or turnover into various analytes, e.g. lactate, based on using α -D-glucose, 6,6- dideutero- α -D-glucose, D_2O , 3-O-methylglucose, 6,6- dideutero- α -D-3-O-trideuteromethylglucose, 6,6- dideutero- α -D-3-O-methylglucose, - α -D-2-O-methylglucose, and other deuterated forms of constituents of nutrient broths and agars associated with the particular cell/tissue types desired. In another embodiment, the method measures indicia of protein metabolism, such as by use of deuterated amino acids.

The method typically comprises contacting the cells of interest with a deuterated agent and probing vibrational spectra emitted by the cells. In one embodiment, the method comprises obtaining a container loaded with deuterated materials, introducing cells into the container whereby the cells are in contact with the deuterated materials, and obtaining vibrational spectra emitted by the cells. The vibrational spectra emitted by the cells are indicative of metabolism, thereby providing an indication of viability of the cells, such that greater metabolic activity is indicative of greater viability. The vibrational spectra can be Raman spectra (including resonance Raman), infrared and/or near

infrared spectra. In some embodiments, the method further comprises subsequently placing the cells in an environment free of the deuterated materials. The rate of clearance of the deuterated material is then monitored to detect viability of the cells. Healthy, viable cells will clear the deuterated material more quickly than nonviable cells.

5

Normalization

The method of the invention can be used for making comparisons between measurements and samples through the use of a "normalizer". In general, any process is more useful with a means for, within a particular sample of cells, controlling or correcting data for the number of cells being probed. A process that relies solely on measuring the intensity of a deuterium related spectral feature by itself is susceptible to errors induced by variation in sample volume, purity and other artifacts of sample preparation and handling. As one example, if at some point the sample of cells is pre-concentrated by centrifugation before spectroscopic probing, then the amount of spectroscopic signal obtained, e.g. Raman or fluorescence or absorption in some spectral region, will depend on how much of the original mother liquor is present or how dispersed the "plug" becomes during sample manipulation/transportation. The amount of the volume of the "plug" that is probed will also depend on the laser or other light source power in a nonlinear and generally complicated manner because there is always a particular geometry of the light delivery and collection optical systems. It is therefore generally useful to use some internal spectroscopic or possibly physical optical signal, i.e. Rayleigh scattered light, that is obtained simultaneously with the deuterated signal, that is indicative of the number of cells being probed or the volume being probed in cases where the density of cells can be independently ascertained.

To cite a few examples of possible normalizers, we note that all living cells contain protein, carbohydrate and lipid in relatively large quantities. Since spontaneous Raman scattering is generally speaking a parts per thousands analytic technique, or equivalently, only capable of providing quantitation down to millimolar concentrations, the choice of normalizer based on a Raman signal, as in the first example given above, is limited to

materials that occur naturally in cells in appropriate concentrations. Normalization can be achieved by determining a ratio between spectral features at a target wavenumber and spectral features at a reference wavenumber. The area under a peak, in spectral counts at or spanning a range of the target wavenumber(s) is divided by the area under a peak in
5 spectral counts at or spanning a range of the reference wavenumber(s). A target or reference wavenumber refers to the spectral feature associated with a particular attribute, such as a C-D bond (at 2400-2600 cm^{-1}).

In the case of vibrational spectroscopies, some spectroscopic features are more appropriate than others for the purpose of volume/cell number normalization. Although
10 the concentration of a given protein may only be micromolar or below, each protein molecule contains thousands to tens of thousands of amide linkages. Each of these bonds produces a spectroscopic signal at nearly the same Raman shift and so the "piling up" of the signal makes the various amide related features, particularly amide I and amide III, strong and often measurable with high signal to noise ratio. Thus the example given
15 above suggested using the so-called amide I signal as a normalizer. Clearly the amide III and other types of amide signals may also be appropriate depending on the specific circumstances of the system in question, e.g. the composition of the nutrient medium and the ability to wash away potentially interfering substances before normalization.

In some cases the normalizer will involve using the CH_2 deformation Raman feature at
20 and around 1450 cm^{-1} . This feature will often be useful because it is a major constituent of phospholipid membranes that are an integral part of any animal cell. Each phospholipid molecule contains around 14-18 of these structural features and so the accumulation of these features allows a convenient and general method for normalization. Note that some nutrient media contain much more protein constituents,
25 i.e. amino acids, than lipid and, in these cases, the adequacy of washing before spectroscopic cell viability assessment is less critical. In the same vein, assessment of plant cell viability will necessarily involve cellulose related features because cell walls are composed of cellulose.

The use of a normalizer is most important when the spectroscopic signals in question are obtained in two separate steps. That is, if the undeuterated spectroscopic signal is obtained using one device and the deuterated signal is obtained using another device, for example because the signals occur in different spectroscopic regions, then a normalizer is appropriate. If the two signals are obtained simultaneously, then a normalizer as described is useful when comparing the results of two separate assessments. But, if the two signals are simultaneously obtained on the same device, then a normalizer is not necessary. Nevertheless, the use of a normalizer may still improve the overall precision of the process.

10 In some embodiments, the 980 nm fluorescence of water is used as a normalizer in concert with Raman excited by a laser having wavelength not shorter than 785 and possibly as long as 980 nm. Choosing the laser wavelength in this range assures that there will be some water fluorescence that can be used to assess water content. The water content can be used as a volume normalizer, i.e. if the cell concentration is known, then a measure of the volume of water in a specific sample can be used as a measure of the number of cells being probed.

In one embodiment of the method, one first incubates the cells in question in an appropriately deuterated medium before removing cells in aliquots to incubate in nondeuterated medium. In this case, the same ratio is measured but now it decreases at a different rate for less viable cells than for more viable cells. In the case of using another analyte in conjunction with glucose (e.g. lactate), there will be more than one Raman feature in the 2200-2500 cm^{-1} range. A ratio is formed involving each of these signals and either the amide I (1600-1700 cm^{-1}) or the amide III (1200-1350 cm^{-1}) as the denominator. The percentage decrease of one feature in the 2200-2500 cm^{-1} range should be matched by an equal percentage increase in the other consistent with the stoichiometry of the overall process.

The amide I Raman feature (1600-1700 cm^{-1}) in the spectra of the cells and tissues and surrounding medium can be used as a normalizer (fluorescence) for the deuterium

specific Raman features to obtain a “specific activity” of the deuterium. The specific activity can be calculated upon probing cell volume, tissue volume or volume of extracellular medium. This quantity will be useful in obtaining optimal quantitative data/analysis for use as a kinetic probe of cell-tissue viability.

- 5 Fluorescence generated by Raman excitation source can be used as a measure of cell volume, tissue volume or volume of extra-cellular medium being probed by Raman features. See, e.g., U.S. Patent No. 6,044,285. In general, porphyrin and long chain terpenoid and similarly derived substances, such as cytochrome, hemoglobin, chlorophyll and the carotenes, respectively, are known to fluoresce in the NIR, thus providing an
10 opportunity to normalize for cell volume without adding exogenous substances, such as fluorophores.

- The amide III Raman feature ($1200\text{--}1350\text{ cm}^{-1}$) in the spectra of the cells and tissues and surrounding medium can be used as a normalizer for the deuterium specific Raman features to obtain a “specific activity” of the deuterium. The specific activity can be
15 calculated upon probing cell volume, tissue volume or volume of extra-cellular medium. This quantity will be useful in obtaining the best quantitative data/analysis for use as a kinetic probe of cell-tissue viability.

- The H_2O fluorescence (980 nm) generated by Raman or a separate excitation source is useful as a measure of cell volume, tissue volume or volume of extra-cellular medium
20 being probed by Raman features.

- The use of the integral of all the Raman features ($300\text{--}1850\text{ cm}^{-1}$) in the spectra of the cells and tissues and surrounding medium can serve as a normalizer for the deuterium specific Raman features to obtain a “specific activity” of the deuterium. The specific activity can be calculated upon probing cell volume, tissue volume or volume of
25 extracellular medium. This quantity will be useful in obtaining optimal quantitative data/analysis for use as a kinetic probe of cell-tissue viability.

Kits

The invention provides kits for carrying out the methods of the invention. In one embodiment, the invention provides a plate with wells, depression(s) or compartment(s) into which a sample of cell or tissue culture can be placed and the spectroscopic process of this invention can be applied. In some embodiments, the plate is preloaded with appropriate deuterated materials to allow multiple kinetic measurements to be made to determine cell/tissue viability. The combination of this process and associated device with the aspects of a plate reader needed to provide cell count information allow the invention to provide absolute calibration to the device. The invention also provides preloaded centrifuge tubes that can be loaded with the cells, processed and then spectroscopically probed in the spun-down tube full of cells. Typically, the centrifuge tubes comprise glass or plastic. In a typical embodiment, the last step in the processing is a spin down step, e.g. following a series of washes.

EXAMPLES

The following examples are presented to illustrate the present invention and to assist one of ordinary skill in making and using the same. The examples are not intended in any way to otherwise limit the scope of the invention.

Example 1: Raman Spectra of Living Cells

This Example demonstrates that the spectroscopic analysis methods can be used to characterize viability of living cells. In this example, phosphate buffered saline (PBS) is used as a control to show that contacting cells with an isosmotic solution containing H₂O (rather than D₂O) does not alter spectroscopic features.

Human leukocytes (HL60-leukemia) were grown by standard tissue culture techniques. An aliquot of the cells was taken from the main culture and exposed to H₂O in phosphate buffered saline (PBS). Raman spectra of the cells were taken at progressively later times to show the effect of the water in the saline solution on the cells. There is clearly very little

effect of adding water, as shown in Figure 1. The largest change observed involved going from the “dry” cells to the water added state. Note, however, that the cells were not actually dry. They had been spun down and there was just a small amount of residual water inside but not on the top of the “plug”. The change in the reflection loss at the top of the cell plug, which varied when the PBS was added, caused the relatively large change in going from “dry” to “wet”.

Example 2: Raman Spectra of Living Cells Using Deuterated Water

This Example demonstrates that the spectroscopic analysis methods can be used to characterize the membrane function of living cells using deuterated agents. Viable cells, with intact membranes, will exhibit features associated with D_2O replacement of H_2O over time. In contrast, nonviable cells will exhibit immediate D_2O replacement. Accordingly, the time course of D_2O exchange can be used as a measure of cell viability.

Figure 2 shows Raman spectra of human leukocytes that were grown by standard tissue culture techniques. The spectra were cut off at a slightly different Raman shift. The right-most peak looks rounded, but the spectra of both “dry” samples (used in Examples 1 and 2) were essentially identical. An aliquot of the cells was taken from the main culture and exposed to D_2O in phosphate buffered saline (PBS). Raman spectra of the cells were taken at progressively later times to show the effect of the deuterated water in the saline solution on the cells. There is clearly an effect. Note that the large change observed going from the “dry” cells (“dry” has same meaning as in Example 1) is the same as above, i.e. the change in the reflection loss at the top of the cell plug varied when either the deuterated or protonated PBS was added.

The large change near 1270 and 960 wavenumbers and other places as well is due to exchange of protons at the amide linkages. As exposure to the D_2O extends from 30 minutes to 60 minutes and then 90 minutes, significant changes are observed as D_2O is exchanged for H_2O . The rate of this exchange depends on the viability of the cells, as cells whose membranes have lost integrity will exhibit immediate changes on the order of changes that would take 90 minutes to occur in healthy, viable cells. The spectra in

Figure 2 clearly show that deuterated agents can be employed in concert with live cells to probe their membrane function.

In order to compare the two sets of spectra quantitatively, normalization is used. For example, one could use the ratio of 1270 cm^{-1} to 1450 cm^{-1} in each set to compare to each other without having to use the 960 cm^{-1} feature. Note also that these spectra do not cover a broad enough range to allow monitoring at $2400\text{-}2600\text{ cm}^{-1}$. Accordingly, one could not simultaneously look at CD bonds and get the spectra shown below unless normalization of some kind is employed. In this case either the 1450 cm^{-1} (CH_2 deformation) or the 1670 cm^{-1} (amide I) could be used for normalization. Note most of the entire feature from $1200\text{ to }1350\text{ cm}^{-1}$ is called amide III. For both sets of spectra, the fluorescence has been subtracted off using a standard technique(101-7 as described in US 6,044,285), resulting in negative counts in some regions. We have produced the same demonstration using Jurket cells, another standard cell line often used in tissue culture.

Those skilled in the art will appreciate that the conceptions and specific embodiments disclosed in the foregoing description may be readily utilized as a basis for modifying or designing other embodiments for carrying out the same purposes of the present invention. Those skilled in the art will also appreciate that such equivalent embodiments do not depart from the spirit and scope of the invention as set forth in the appended claims.